

CARDIOVASCULAR CANNULATION PRODUCTS	3-9
Arterial cannulae	3
Aortic catheter	7
Flex line venous cannulae	8
Flex line two stage venous cannulae	9
CARDIOPLEGIA PRODUCTS	10 - 16
Aortic root cannulae	10
Retrograde cardioplegia cannula	12
Cardioplegia needle	12
Cardioplegia adapters	13
Ostial perfusion cannulae	16
VENT CATHETERS	17 - 20
Silicone vent catheter	17
PVC vent catheters	18
SUCKERS	21-23
Rigid suckers	21
Flexible suckers	21
Yankauer suction tube	22
Pericardial sump	23
GAS DIFFUSER	24
HEMOCONCENTRATORS	25
ACCESSORIES	26 - 29
Suction connecting tubes	26
Aspiration tubes	26
Porous vent plugs	26
Vented connector cap	26
Vessel cannulae	27
Tubing organizer	27
Connectors	28
Cannulation tourniquet sets	29
Quick prime line	29
Pressure monitoring line	29
MAMMARIA RETRACTOR	30





REINFORCED, CURVED TIP WITH FLANGE

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP
	20 Fr	33cm	A2016	A2016V		24 cm	B2016	B2016V	
WITHOUT CONNECTOR	22 Fr	33cm	A2216	A2216V		24 cm	B2216	B2216V	
	24 Fr	33cm	A2416	A2416V		24 cm	B2416	B2416V	
	20 Fr	33cm	A20161		A20161V	24 cm	B20161		B20161V
WITH CONNECTOR	22 Fr	33cm	A22161		A22161V	24 cm	B22161		B22161V
	24 Fr	33cm	A24161		A24161V	24 cm	B24161		B24161V
	20 Fr	33cm	A20162		A20162V	24 cm	B20162		B20162V
WITH LUER CONNECTOR	22 Fr	33cm	A22162		A22162V	24 cm	B22162		B22162V
	24 Fr	33cm	A24162		A24162V	24 cm	B24162		B24162V





REINFORCED, STRAIGHT TIP WITH RING

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- √ Kink resistant reinforced tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP
	20 Fr	33 cm	A2013	A2013V		24 cm	B2013	B2013V	
WITHOUT CONNECTOR	22 Fr	33 cm	A2213	A2213V		24 cm	B2213	B2213V	
	24 Fr	33 cm	A2413	A2413V		24 cm	B2413	B2413V	
	20 Fr	33 cm	A20131		A20131V	24 cm	B20131		B20131V
WITH CONNECTOR	22 Fr	33 cm	A22131		A22131V	24 cm	B22131		B22131V
	24 Fr	33 cm	A24131		A24131V	24 cm	B24131		B24131V
	20 Fr	33 cm	A20132		A20132V	24 cm	B20132		B20132V
WITH LUER CONNECTOR	22 Fr	33 cm	A22132		A22132V	24 cm	B22132		B22132V
	24 Fr	33 cm	A24132		A24132V	24 cm	B24132		B24132V





NON-REINFORCED, CURVED TIP WITH FLANGE

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- √ Kink resistant pvc tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
	20 Fr	25 cm	A2006	
WITHOUT CONNECTOR	22 Fr	25 cm	A2206	
COMMECTOR	24 Fr	25 cm	A2406	
	20 Fr	25 cm	A20061	A20061V
WITH CONNECTOR	22 Fr	25 cm	A22061	A22061V
	24 Fr	25 cm	A24061	A24061V
	20 Fr	25 cm	A20062	A20062V
WITH LUER CONNECTOR	22 Fr	25 cm	A22062	A22062V
	24 Fr	25 cm	A24062	A24062V





NON-REINFORCED, STRAIGHT TIP WITH RING

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant pvc tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
	20 Fr	25 cm	A2003	
WITHOUT CONNECTOR	22 Fr	25 cm	A2203	
	24 Fr	25 cm	A2403	
	20 Fr	25 cm	A20031	A20031V
WITH CONNECTOR	22 Fr	25 cm	A22031	A22031V
	24 Fr	25 cm	A24031	A24031V
	20 Fr	25 cm	A20032	A20032V
WITH LUER CONNECTOR	22 Fr	25 cm	A22032	A22032V
	24 Fr	25 cm	A24032	A24032V





AORTIC CATHETER

Andocor aortic catheters are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth beveled tip
- ✓ Movable suture ring



	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
	20 Fr	24 cm	AC201341	AC201341V
	22 Fr	24 cm	AC221341	AC221341V
WITH CONNECTOR	24 Fr	24 cm	AC241341	AC241341V
	26 Fr	24 cm	AC261341	AC261341V
	30 Fr	24 cm	AC301341	AC301341V
	22 Fr	24 cm	AC221342	AC221342V
WITH LUER	24 Fr	24 cm	AC241342	AC241342V
CONNECTOR	26 Fr	24 cm	AC261342	AC261342V
	30 Fr	24 cm	AC301342	AC301342V
SPECIAL CONFIGURATIO	N: HARD	TIP		
WITH CONNECTOR	24 Fr	24 cm	AC241341H	AC241341VH
WITH LUER CONNECTOR	24 Fr	24 cm	AC241342H	AC241342VH





REINFORCED WITH LIGHTHOUSE TIP

FLEX LINE VENOUS CANNULAE


Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.





Features & Benefits

- ✓ Single piece construction
- ✓ Open and flexible lighthouse tip
- ✓ Easy positioning
- ✓ Insertion depth marker



SIZE	STRAIGHT TIP	RIGHT ANGLED 90° TIP	CONNECTION SITE
20 Fr	01V201L8	01V201L9	3/8"
22 Fr	01V221L8	01V221L9	3/8"
24 Fr	01V241L8	01V241L9	3/8"
26 Fr	01V261L8	01V261L9	3/8"
28 Fr	01V281L8	01V281L9	3/8"
30 Fr	01V301L8	01V301L9	3/8"
32 Fr	01V321L8	01V321L9	3/8"
34 Fr	01V341L8	01V341L9	3/8"
36 Fr	01V361L8	01V361L9	1/2"

SIZE	NON-REINFORCED	CONNECTION SITE
32 Fr	VC3208	3/8"





PROXIMAL AND DISTAL REINFORCED WITH LIGHTHOUSE TIP

FLEX LINE TWO STAGE VENOUS CANNULAE

Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.




Features & Benefits

- ✓ Single piece construction
- ✓ Open and flexible lighthouse tip
- ✓ Easy positioning
- ✓ Obturator
- \checkmark Insertion depth marker



DETAILS	
Connection site	1/2"
Length	40 cm
Qty/box	20 Pcs



	SIZE	STANDARD BODY	FLAT BODY
	32/40 Fr	01V32L40	01V32L40FB
WITHOUT	34/46 Fr	01V34L46	01V34L46FB
CONNECTOR	36/46 Fr	01V36L46	01V36L46FB
	36/51 Fr	01V36L51	01V36L51FB
	32/40 Fr	01V32L401	01V32L401FB
WITH	34/46 Fr	01V34L461	01V34L461FB
CONNECTOR	36/46 Fr	01V36L461	01V36L461FB
	36/51 Fr	01V36L511	01V36L511FB
WITH LUER	32/40 Fr	01V32L402	01V32L402FB
	34/46 Fr	01V34L462	01V34L462FB
CONNECTOR	36/46 Fr	01V36L462	01V36L462FB
	36/51 Fr	01V36L512	01V36L512FB








DETAILS	
Length	15 cm
Qty/box	20 Pcs

WITH OR WITHOUT VENT LINE

AORTIC ROOT CANNULAE

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart.

The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.



	CODE	TIP SIZE	NEEDLE SIZE
WITHOUT VENT LINE	ARC07	6 Fr	16 G
	ARC04	7 Fr	14 G
	ARC01	9 Fr	12 G
WITH VENT LINE	ARC08	6 Fr	16 G
	ARC05	7 Fr	14 G
	ARC02	9 Fr	12 G



WITH ANTI-BLEEDING VALVE



Features & Benefits





	CODE	TIP SIZE	NEEDLE SIZE
WITHOUT VENT LINE	ARC07V	6 Fr	16 G
	ARC04V	7 Fr	14 G
	ARC01V	9 Fr	12 G
WITH VENT LINE	ARC08V	6 Fr	16 G
	ARC05V	7 Fr	14 G
	ARC02V	9 Fr	12 G

SPECIAL CONFIGURATIONS





	CODE	TIP SIZE	NEEDLE SIZE	LENGTH
LONG TYPE WITHOUT VENT LINE	ARC14	7 Fr	14 G	28 cm
WITH VENT LINE WITH 2 CLAMPS	ARC10	7 Fr	14 G	15 cm






Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedures.



Features & Benefits

- ✓ Self-inflating
- ✓ Silicone balloon
- ✓ Malleable stylet



DETAILS		
Code	RCC01	
Length	31 cm	
Balloon diameter	18 mm	
Size	14 Fr	
Qty/box	20 Pcs	

CARDIOPLEGIA NEEDLE

Andocor cardioplegia needles are intended for use in delivery of cardioplegia solution during cardiopulmonary bypass surgery.



Features & Benefits

- ✓ Protective cap
- ✓ Stainless steel needle
- ✓ Luer lock connection



DETAILS		
Code	CN16	
Length needle	1,4 cm	
Size needle	16G	
Total length	5 cm	
Qty/box	20 Pcs	



CARDIOPLEGIA ADAPTERS

Andocor cardioplegia adapters connect to the aortic root cannulae or vessel cannulae and are intended for the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass.

	DETAILS	
Qty/box		20 Pcs

DIFFERENT CONFIGURATIONS				
CSY03	PERFUSION ADAPTER WITH 2 MALE LUER LOCK AND 1 FEMALE LUER LOCK CON- NECTOR AND COLOR- CODED CLAMPS			
CSY04	MULTIPLE PERFUSION ADAPTER FOR AORTIC ROOT AND THREE VEIN GRAFTS			
CSY05	MULTIPLE PERFUSION ADAPTER FOR AORTIC ROOT AND THREE VEIN GRAFTS WITH ATTACHED VESSEL CANNULAE			
CSY06	DUAL PERFUSION SET WITH MALE LUER LOCK, CONE CONNECTOR AND COLOR CODED CLAMPS	6		
CSY07	MULTIPLE PERFUSION ADAPTER WITH VENT LINE AND ATTACHED VESSEL CANNULAE			



CSY08	AORTIC PERFUSION SET WITH VENT LINE, TOTAL LENGTH = 30 CM	•
CSY09	DUAL PERFUSION SET WITH MALE LUER LOCK, FEMALE LUER LOCK CONNECTORS AND COLOR-CODED CLAMPS	
CSY10	MULTIPLE SET FOR PERFUSION OF THREE VEIN GRAFTS	***
CSY11	MULTIPLE SET FOR PERFUSION OF THREE VEIN GRAFTS WITH ATTACHED VESSEL CANNULAE	**
CSY12	MULTIPLE PERFUSION ADAPTER FOR RETROGRADE AND AORTIC ROOT WITH VENT LINE	
CSY13	AORTIC PERFUSION SET WITH VENT LINE, TOTAL LENGTH = 18 CM	6
CSY14	AORTIC PERFUSION SET WITH 1/4"x1/4"x1/4" CONNECTOR	
CSY15	STRAIGHT ADAPTER. 14 CM	



CSY16	STRAIGHT ADAPTER. 15 CM	
CSY17	PERFUSION ADAPTER 1/4"	
CSY18	PERFUSION ADAPTER 3/8"	-ME
CSY19	PERFUSION ADAPTER 1/4" WITH 3-WAY STOPCOCK	
CSY20	PERFUSION ADAPTER WITH 2 FEMALE LUER LOCKS, 1 MALE LUER LOCK AND COLOR-CODED CLAMPS	
CSY21	RECIRCULATION Y-ADAPTER	
CSY22	MULTIPLE PERUSION ADAPTER WITH VENT LINE	
CSY23	SIX MULTIPLE PERFUSION ADAPTER	



BALLOON TIP



TIP SIZE

3 mm

3 mm

DETAILS			
Length	30 cm		
Qty/box	20 Pcs		

(

Features & Benefits

- ✓ Self-inflating silicone balloon
- √ Female luer lock connection

BASKET TIP

CODE

OC06A

OC07A

OSTIAL PERFUSION CANNULAE

BALLOON SIZE

6 mm

7 mm

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

DETAILS		
Length	21 cm	
Qty/box	20 Pcs	

Features & Benefits

- ✓ Basket tip
- ✓ Different configurations
- ✓ Malleable stainless steel shaft

	TIP SIZE	CODE	
	10 Fr	OC1045	
WITH 45° TIP	12 Fr	OC1245	
	14 Fr	OC1445	
	10 Fr	OC1090	
WITH 90° TIP	12 Fr	OC1290	
	14 Fr	OC1490	

VENT CATHETERS





SILICONE VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



Features & Benefits

- ✓ With guidewire stylet
- ✓ Soft and flexible silicone tubing
- ✓ Soft tip



DETAILS				
Code	SV20			
Size	20 Fr			
Length	35 cm			
Tip perforation	6,5 cm			
Connection site	1/4"			
Qty/box	20 Pcs			

DETAILS				
Code	SV16			
Size	16 Fr			
Length	34 cm			
Tip perforation	6,5 cm			
Connection site	1/4"			
Qty/box	20 Pcs			





PVC VENT CATHETER

Andocor vent catheters are medical devices used to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.

DETAILS			
Connection site	1/4"		
Qty/box	20 Pcs		



Features & Benefits

- √ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker

	SIZE	TIP PERFORATION	LENGTH	CODE
WITHOUT CTVLET	13 Fr	3,8 cm	32 cm	PVBA13
WITHOUT STYLET	20 Fr	10 cm	36,5 cm	PVBA20

VENT CATHETERS



OPEN TIP, STRAIGHT



	SIZE	TIP PERFORATION	LENGTH	CODE
WITH PLASTIC STYLET	17 Fr	2,5 cm	32 cm	PVSOR170
	18 Fr	2,5 cm	32 cm	PVSOR180
	17 Fr	6 cm	32 cm	PVSOR171
	18 Fr	6 cm	32 cm	PVSOR181
	17 Fr	9 cm	32 cm	PVSOR172
	18 Fr	9 cm	32 cm	PVSOR182




	SIZE	TIP PERFORATION	LENGTH	CODE
WITH GUIDEWIRE STYLET	18 Fr	6 cm	32 cm	PVSOG181





CLOSED TIP, STRAIGHT



	SIZE	TIP PERFORATION	LENGTH	CODE
	17 Fr	2,5 cm	32 cm	PVSBG170
WITH GUIDEWIRE STYLET	18 Fr	2,5 cm	32 cm	PVSBG180
	17 Fr	6 cm	32 cm	PVSBG171
	18 Fr	6 cm	32 cm	PVSBG181
	20 Fr	6 cm	32 cm	PVSBG201
	17 Fr	9 cm	32 cm	PVSBG172
	18 Fr	9 cm	32 cm	PVSBG182





	SIZE	TIP PERFORATION	LENGTH	CODE
	17 Fr	2,5 cm	32 cm	PVSBR170
WITH PLASTIC STYLET	18 Fr	2,5 cm	32 cm	PVSBR180
	17 Fr	6 cm	32 cm	PVSBR171
	18 Fr	6 cm	32 cm	PVSBR181
	20 Fr	6 cm	32 cm	PVSBR201
	17 Fr	9 cm	32 cm	PVSBR172
	18 Fr	9 cm	32 cm	PVSBR182





Andocor suckers are used for the removal of excess fluid from the surgical field.

DETAILS			
Code	RS01		
Length	33 cm		
Qty/box	20 Pcs		



Features & Benefits

- ✓ Metal suction tube with plastic handle
- ✓ Easy & comfortable directing



FLEXIBLE & MALLEABLE SUCKERS

Andocor suckers are used for the removal of excess fluid from the surgical field.

DETAILS		
Connection site	1/4"	
Qty/box	20 Pcs	



Features & Benefits

- ✓ Provide continuous suction
- ✓ Flexible or malleable tubing

TUBING	CODE	LENGTH	FRENCH SIZE
FLEXIBLE	FS300	36 cm	17 Fr
MALLEABLE	FS301	45 cm	17 Fr





YANKAUER SUCTION TUBE

DETAILS

Qty/box 20 Pcs

Andocor yankauer suction tubes are intended to aspirate body waste fluid during surgery.

CONFIGURATONS	TIP	CODE
VENTED	OPEN TIP	ST04
VENTED	BULB TIP	ST02
NON-VENTED	OPEN TIP	ST03
	BULB TIP	ST01





Andocor suckers are used for the removal of excess fluid from the surgical field.



Features & Benefits

- ✓ Smooth and flush design
- ✓ Weighted tip
- ✓ Drainage of pericardial sack





METAL TIP

PS40M | PS40ML | PS40MLH

PS40M1

CODE	DESCRIPTION	TUBING LENGTH	EXTERNAL DIAMETER TIP
PS40M	METAL TIP	36 cm	6 mm
PS40ML	LONG TUBING	80 cm	6 mm
PS40MLH	LONG, HARD TUBING	80 cm	6 mm
CODE	DESCRIPTION	TUBING LENGTH	EXTERNAL DIAMETER TIP
PS40M1	LARGE METAL TIP	36 cm	7,8 mm

WIREWOUND TIP



DETAILS			
Code	PS40		
Size	16 Fr		
Length	36 cm		
Qty/box	20 Pcs		





GAS DIFFUSER

Andocor gas diffuser is a medical device for intraoperative wound ventilation with co_2 to prevent air embolism and bacteriological contamination during surgery. By creating a continuous flow of co_2 in the intraoperative wound, air and bacteriological particles are removed.



Features & Benefits

- ✓ Innovative microporous diffuser tip
- ✓ Hydrophobic diffuser tip remains functional even when wet
- ✓ Bacterial filter
- ✓ Malleable section allows precise placement



WITH MALLEABLE TUBING

DETAILS			
Code	GAS225		
Malleable section	40 cm		
Length	225 cm		
Qty/box	20 Pcs		

SPECIAL CONFIGURATION

DOUBLE PACKED GAS DIFFUSER

DETAILS			
Code	GAS225DP		
Malleable section	40 cm		
Length	225 cm		
Qty/box	20 Pcs		



WITH FLEXIBLE TUBING

DETAILS			
Code	GAS225F		
Flexible section	30 cm		
Length	225 cm		
Qty/box	20 Pcs		

HEMOCONCENTRATORS





HEMOCONCENTRATORS

Andocor hemoconcentration filters are intended to relieve or mitigate overhydratation in patients undergoing cardiopulmonary procedures and to increase concentrations of cells and proteins in the blood.



Features & Benefits

- ✓ Pre-rinsed or non-rinsed filters
- ✓ Standard or custom-made blood lines

DESCRIPTION	FILTERS	CODE	QTY/BOX
HEMOCONCENTRATOR	NON-RINSED	03HC0250	10 Pcs
HC0250	PRE-RINSED	03HP0250	
HEMOCONCENTRATOR	NON-RINSED	03HC0600	6 Pcs
HC0600	PRE-RINSED	03HP0600	
HEMOCONCENTRATOR	NON-RINSED	03HC1000	6 Pcs
HC1000	PRE-RINSED	03HP1000	
TUBING SET		03HCS001	40 Pcs

TECHNICAL FEATURES			
	HC0250	HC0600	HC1000
INTENDED USE	PAEDIATRIC	ADULT	ADULT
MEMBRANE AREA	0.25 m ²	0.6 m ²	1.0 m ²
PRIMING BLOOD VOLUME	18 ml	51 ml	76 ml
MAX. TRANS- MEMBRANE PRESSURE	500 mmHg	500 mmHg	500 mmHg
UNIT LENGTH	180 mm	180 mm	265 mm
ULTRAFILTRATION RATE AVERAGE KUF IN VIVO 22 ml/h/mmHg/m²	5.5	13.2	22.0



SUCTION CONNECTING TUBES

Andocor suction connecting tubes are used for removal of excess fluid from a surgical field.



Features & Benefits

- ✓ Double female funnel
- ✓ Double packed

CODE	SIZE	LENGTH	QTY/BOX
TUB69180	6x9 mm	180 cm	20 Pcs
TUB69380	6x9 mm	380 cm	10 Pcs



ASPIRATION TUBES

Andocor aspiration tubes are used as an accessory for blood aspiration during cardiovascular surgery and cannulation.

CODE	LENGTH	QTY/BOX
ASP180	180 cm	25 Pcs
ASP380	380 cm	25 Pcs



POROUS VENT PLUG

CODE	SIZE	QTY/BOX
VP38	3/8"	40 Pcs



VENTED CONNECTOR CAP

CODE	SIZE	QTY/BOX
VCC01	3/8"	40 Pcs



ACCESSORIES









VESSEL CANNULAE

Andocor vessel cannulae are used to help check for leaks in a harvested vein that will be used in a graft with a sterile saline solution.

Additionally, the vessel cannula can be used for antegrade cardioplegia administration.

	CODE	LENGTH	QTY/BOX
TRANSPARANT, FLEXIBLE	VC3	5 cm	40 Pcs
WITH CLOSED BACK CHECK VALVE	VC3B	6,5 cm	40 Pcs
WHITE, RIGID	VC04	5 cm	40 Pcs



TUBING ORGANIZER

DETAILS			
CODE	TO10		
SLOTS & TYPE	1 TUBE 1/8"		
	6 TUBES 1/4"		
	2 TUBES 3/8"		
	1 TUBE 1/2"		
QTY/BOX	10 Pcs		







CONNECTORS

DETAILS			
Qty/box	20 pcs		

Andocor connectors are accessories for the connection of cannulae to the extracorporeal cardiopulmonary bypass system and to connect tubing.


		CODE	OUTLET 1	OUTLET 2	OUTLET 3
		04CS001	1/2"	1/2"	
		04CS002	3/8"	3/8"	
		04CS003	1/4"	1/4"	
	WITHOUT LUER	04CS004	3/8"	1/2"	
	LOCK	04CS005	1/4"	3/8"	
		04CS006	3/16"	1/4"	
CTD ALCUT		04CS007	3/16"	3/16"	
STRAIGHT CONNECTOR		04CS008	3/8"	3/16"	
CONNECTOR		04CS0012	1/2"	1/2"	
		04CS0022	3/8"	3/8"	
	WITH	04CS0032	1/4"	1/4"	
	LUER LOCK	04CS0042	3/8"	1/2"	
		04CS0052	1/4"	3/8"	
		04CS0062	3/16"	1/4"	
		04CS0072	3/16"	3/16"	
		04CY001	3/8"	3/8"	3/8"
		04CY002	1/4"	1/4"	1/4"
	\\/\!\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\	04CY003	3/8"	3/8"	1/2"
	WITHOUT LUER LOCK	04CY004	1/4"	1/4"	3/8"
Y-CONNECTOR	LOCK	04CY005	1/2"	1/2"	1/2"
		04CY006	3/8"	1/4"	3/8"
		04CY007	1/2"	1/4"	1/2"
	WITH	04CY0012	3/8"	3/8"	3/8"
	LUER LOCK	04CY0022	1/4"	1/4"	1/4"



CANNULATION TOURNIQUET SETS

Andocor cannulation tourniquet sets are used to secure arterial and venous cannulae during cardiopulmonary bypass.


CODE	SNARE	RED TUBE	BLUE TUBE	CLEAR TUBE	QTY/BOX
CT01	1	1	1	3	20 Pcs
CT02	1	1	2	3	20 Pcs
CT03	2	1	2	5	20 Pcs
CT04	1	1	1		20 Pcs
CT05	1	2	1		20 Pcs
CT06	2	2	2	4	20 Pcs
CT07	2	4	2	2	20 Pcs
CT08	1	5			20 Pcs
CT09	1			3	20 Pcs



QUICK PRIME LINE

Andocor quick prime lines facilitate the transfer of fluids during the priming of the extracorporeal circuit.

DETAILS			
Code	QPL100		
Length	100 cm		
Qty/box	20 pcs		



PRESSURE MONITORING LINE

Andocor pressure monitoring lines are used for extending medical tubing and pressure monitoring.

CODE	TUBE END	LENGTH	QTY/BOX	PRESSURE
PML180MM	MALE/MALE	180 cm	20 Pcs	800 psi
PML180MF	MALE/FEMALE	180 cm	20 Pcs	800 psi








MAMMARIA RETRACTOR

The mammaria retractor is used for the preparation of the internal mammary artery and can be utilised for patients weighing over 25 kg (55 lbs).

DETAILS

SURGICAL STAINLESS STEEL

ASYMMETRICAL OPENING UP TO 70 mm MAX



ANDOCOR n.v.

Kwikaard 104

B-2980 Zoersel

Belgium

T+32 (0)3 380 15 95

F +32 (0)3 380 15 90

info@andocor.com

www.andocor.com